

molecular biology

applied biosystems

invitrogen



Prepare for discovery

Molecular biology workflow solutions

Find it at fishersci.com

 **fisher** scientific
part of Thermo Fisher Scientific



AmpliFly

AmpliFly



Prepare for discovery with molecular biology

In molecular biology research, you may take different routes to reach your destinations. Choices can often be confusing, and the path you take may impact the experience and overall success of your experiments.

This handbook is intended to guide you by providing technical information and clear choices across the molecular biology workflow to help you soar in your research. Applied Biosystems™ and Invitrogen™ products incorporate the latest innovations to enable quicker results, more assurance, and less optimization in your lab studies.

Prepare to step into discovery, from sample preparation to reverse transcription, PCR, and cloning. Explore this handbook for more direct routes and a first-class experience in your research from sample to result.

Find additional information at thermofisher.com/amplify

Contents



Sample preparation

| | |
|---------------------|----|
| Methods overview | 7 |
| DNA isolation | 9 |
| RNA isolation | 11 |
| Automation platform | 12 |



Reverse transcription

| | |
|---------------------|----|
| Considerations | 15 |
| Reagent selection | 16 |
| Genomic DNA removal | 17 |
| Primers | 18 |



PCR

| | |
|-------------------------------|----|
| Thermal cycler considerations | 21 |
| Plastics essentials | 24 |
| Enzyme properties | 26 |
| Oligo design and selection | 30 |



Electrophoresis

| | |
|--------------------------|----|
| Workflow | 33 |
| E-Gel selection | 34 |
| Electrophoresis reagents | 36 |



Cloning

| | |
|----------------------------|----|
| Technologies overview | 39 |
| Restriction enzyme cloning | 40 |
| PCR cloning | 42 |
| Cloning with synthetic DNA | 44 |
| Transformation | 45 |




Resources

| | |
|--------------------------|----|
| Educational resources | 49 |
| Mobile apps | 50 |
| Custom and OEM solutions | 50 |
| FAQs | 51 |
| Ordering information | 53 |




Sample preparation



Nucleic acid isolation is a crucial first step in the molecular biology workflow, whether you are isolating genomic DNA (gDNA) or RNA. Selecting nucleic acid purification products that are optimized to provide maximum yield, purity, and integrity from virtually any sample type and application is important for your research success.

Find technical resources on nucleic acid isolation at
[thermofisher.com/prepon](https://www.thermofisher.com/prepon)

Benefits and underlying principles of common nucleic acid isolation methods



Organic extraction: Phenol-chloroform solution (e.g., Invitrogen™ DNAzol™ and TRIzol™ reagents)

After homogenizing the sample with TRIzol Reagent, chloroform is added, and the mixture separates into a clear upper aqueous layer containing RNA, an interphase layer, and a pink lower organic layer containing the DNA and protein. RNA is precipitated from the upper aqueous layer with isopropanol. DNA is precipitated from the interphase and organic layers with ethanol. Protein is precipitated from the phenol-ethanol supernatant with isopropanol.

Benefits:

- Efficient lysis of cells and tissue
- Rapid denaturation of nucleases
- Stabilization of nucleic acids
- Great for fatty and cartilaginous samples




Spin columns: Glass fiber, derivatized silica, or ion exchange membrane in column (e.g., Thermo Scientific™ GeneJET™ and Invitrogen™ PureLink™ kits)

Samples are lysed and passed through the membrane using centrifugal or vacuum force. Wash and elution solutions are subsequently passed through the membrane, and the sample is collected into a tube by centrifugation.

Benefits:

- Convenience
- Ease of use
- Throughput flexibility
- Specialized equipment not required



Magnetic beads: 0.5–1.0 µm particles with a paramagnetic core and modified shell (e.g., Applied Biosystems™ MagMAX™ kits and Invitrogen™ Dynabeads™ magnetic beads)


Samples are lysed in solution and allowed to bind nucleic acid to magnetic particles based on specific surface modifications. Application of an external magnetic field rapidly collects the particles. Rounds of release, washes, and recapture enable purification of the desired nucleic acid.

Benefits:

- No risk of clogging
- Increased target capture efficiency
- Rapid collection and concentration of sample
- Specialized equipment not required
- Scalability

Find out more at thermofisher.com/sampleprep

Tools for success in nucleic acid isolation



DNA type

For gDNA extraction, plasmid isolation, and DNA cleanup



RNA type

For purification of total RNA, transcriptome RNA, messenger RNA (mRNA), microRNA (miRNA) and other small RNA, and sequence-specific RNA capture

Find out more at thermofisher.com/prepon

Selecting the right DNA isolation kits

To suit your specific DNA isolation needs, a comprehensive portfolio is available offering three highly developed purification technologies: silica membrane, anion-exchange resin, and switchable surface charge. Refer to the following tables to determine which purification kit is appropriate for your DNA type and sample type. To use our online kit selection guide, go to thermofisher.com/dnaselection

Applied Biosystems™ and Invitrogen™ technologies for gDNA isolation

| Throughput/format | Low throughput or manual (organic) | Medium throughput and spin-column technology | High throughput and 96-well filter plate | High throughput and magnetic bead technology |
|--|------------------------------------|--|--|--|
| Kits | DNAzol reagents | PureLink kits | PureLink Pro gDNA kits | MagMAX DNA kits and Dynabeads products |
| Tissue | DNAzol Reagent | PureLink Genomic DNA Mini Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX DNA Multi-Sample Kit |
| Cells | DNAzol Reagent | PureLink Genomic DNA Mini Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX DNA Multi-Sample Kit |
| Blood | DNAzol BD Reagent | PureLink Genomic DNA Mini Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX DNA Multi-Sample Ultra 2.0 Kit |
| Plant | Plant DNAzol Reagent | PureLink Genomic Plant DNA Purification Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX Plant DNA Isolation Kit |
| Buccal swabs | Not recommended | PureLink Genomic DNA Mini Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX DNA Multi-Sample Ultra 2.0 Kit |
| Bacteria | DNAzol Reagent | PureLink Microbiome DNA Purification Kit | PureLink Pro 96 Genomic DNA Purification Kit | MagMAX DNA Multi-Sample Ultra 2.0 Kit, Dynabeads MyOne Streptavidin products |
| Virus | Not recommended | PureLink Viral RNA/DNA Mini Kit | PureLink Pro 96 Viral RNA/DNA Purification Kit | MagMAX Viral Nucleic Acid Kit, Dynabeads MyOne Streptavidin products |
| Compatibility | | | | |
| Scalable and automatable | | | • | • |
| Thermo Scientific™ KingFisher™ instrument | | | | • |
| qPCR | • | • | • | • |
| NGS | • | • | • | • |

Find out more at thermofisher.com/gdnaprep

Selecting the right DNA purification kits (cont.)

Comparison of DNA cleanup solutions

| DNA cleanup application | PCR cleanup | PCR cleanup | PCR cleanup | Gel extraction | PCR cleanup and gel extraction | PCR cleanup | Sequencing reaction cleanup |
|-------------------------|---|--|---|---|--|---|---|
| Kit | PureLink PCR Purification Kit | PureLink Pro 96 PCR Purification Kit | PureLink PCR Micro Kit | PureLink Quick Gel Extraction Kit | PureLink Quick Gel Extraction Kit and PCR Purification Combo Kit | ChargeSwitch-Pro PCR Clean-Up Kit | Centri-Sep Spin Columns |
| Format | Silica spin/vacuum column | 96-well silica plate | Silica spin column | Silica spin/vacuum column | Silica spin/vacuum column | Derivatized spin/vacuum column | Spin column |
| Product size | <ul style="list-style-type: none"> • 50 preps • 250 preps | <ul style="list-style-type: none"> • 4 plates (4 x 96 rxns) | <ul style="list-style-type: none"> • 10 preps • 50 preps • 250 preps | <ul style="list-style-type: none"> • 50 preps • 250 preps | <ul style="list-style-type: none"> • 50 preps | <ul style="list-style-type: none"> • 10 preps • 50 preps • 250 preps | <ul style="list-style-type: none"> • 100 columns • 32 columns |
| Time | <15 | 20 | ≤10 | <30 | 10–30 | <10 | <5 |
| Elution volume | 50 | 50–150 | 5–20 | 30–100 | 30–100 | 50 | 20 |
| Recovery | >80% | >70% | >80% | Up to 95% | Gel cleanup: >80% PCR cleanup: >95% | NA | NA |
| Primer removal | >99% | NA | >95% | NA | >99% | NA | >98% |

Thermo Scientific™ and Invitrogen™ technologies for plasmid DNA isolation

| Purity grade | Molecular | Transfection | Transfection | Advanced transfection |
|------------------------|---|--|---|--|
| Kits | GeneJET kits | PureLink HiPure kits | PureLink Fast Low-Endotoxin kits | PureLink Expi Endotoxin-Free kits |
| Endotoxin level | Standard (>10 EU/μg) | Low endotoxin (1–10 EU/μg) | Low endotoxin (0.1–1 EU/μg) | Endotoxin-free (<0.1 EU/μg) |
| Yield | 20 μg–1 mg | 20 μg–15 mg | 0.4 mg (midi), 1.5 mg (maxi) | 1.5–15 mg |
| Technology | Silica membrane | Anion exchange (resin) | Advanced silica membrane | Anion exchange (membrane) |
| Total protocol time | 15–60 min | 30–120 min | 30 min | 90–120 min |
| Prep size | Mini, midi, and maxi | Mini, midi, maxi, and giga | Midi and maxi | Maxi, mega, and giga |
| Downstream application | <ul style="list-style-type: none"> • PCR • Nucleic acid labeling • Cloning (digestion, ligation) • Sequencing | <ul style="list-style-type: none"> • Standard transfection • All molecular biology applications • <i>In vitro</i> transcription | <ul style="list-style-type: none"> • Standard transfection • Transfection of certain sensitive cell lines | <ul style="list-style-type: none"> • Primary and stem cell transfection • Gene therapy and vaccine (<i>in vivo</i>) research • Microinjection • All molecular biology applications |

Selecting the right RNA isolation kits

For the quality and performance you need, a full suite of products for RNA isolation is available for a wide range of sample types, throughputs, and input quantities. To use our online kit selection guide, go to thermofisher.com/rnaselection

Applied Biosystems™ and Invitrogen™ technologies for total RNA isolation

| Capabilities | Process a large amount of tissue | Fast isolation of RNA from a variety of samples | High-throughput purification of RNA and DNA | Process cells for gene expression |
|----------------------------|---|---|---|-----------------------------------|
| Kits | TRIzol reagents | PureLink kits | MagMAX kits | Cells-to-C _T kits |
| Prep time | 30–60 min | <20 min | 45 min | ≤10 min |
| Sample types | Most samples, particularly those more difficult to lyse | Bacteria, liquid, blood, cells, yeast, plants, tissue | Cells, blood, plants | Cultured cells |
| Starting material | 100 mg of tissue or 10 ⁷ cells | 10 ⁸ cells, 200 mg of tissue, 250 mg of plant tissue, 0.2 mL of blood, 5 x 10 ⁸ yeast, 10 ⁹ bacteria | 100 mg of tissue or 5 x 10 ⁶ cells | 1–100,000 cells |
| Yield | 10 ⁶ epithelial cells: 8–15 µg 100 mg tobacco leaf: 73 µg (Variable depending on sample) | Up to 350 µg | Variable depending on sample | NA |
| High-throughput compatible | | Yes | Yes | Yes |
| Technology | Organic extraction | Silica membrane spin column/filter plate | Magnetic beads | Crude lysate |



Helpful tip

If you are not ready to process your RNA sample, simply store it in Invitrogen™ RNA/*later*™ Stabilization Solution for use at a later time. Visit thermofisher.com/stabilizerna



Find out more at thermofisher.com/rnapreps

Automation platform

Optimize and automate your DNA and RNA purification workflow with Thermo Scientific™ KingFisher™ purification systems. When used with compatible bead-based reagents, such as MagMAX kits and Dynabeads products, these instruments enable versatile automation of DNA and RNA isolation procedures. Learn more and request a demo at thermofisher.com/kingfisherdemo



Thermo Scientific™ KingFisher™ Duo Prime Purification System

Automated nucleic acid purification of 6 or 12 samples using magnetic bead technology.

Thermo Scientific™ KingFisher™ Flex Purification System

Highly versatile and automated nucleic acid purification of 24 or 96 samples per run using magnetic bead technology.



Thermo Scientific™ KingFisher™ Presto Purification System

Automated nucleic acid purification of 24 to 96 samples per run using magnetic bead technology. This instrument is specifically designed to be paired with a liquid handler for high-throughput, fully automated nucleic acid purification.



Resources

Navigate through the DNA and RNA support categories below to obtain relevant technical information, view tips and tricks when starting an experiment, and find answers to everyday problems.

thermofisher.com/napsupport

thermofisher.com/technicalresources

thermofisher.com/prepforsuccess

thermofisher.com/rnabasics

thermofisher.com/rnahandlingtips

thermofisher.com/magmax

thermofisher.com/dynabeads

Find out more at thermofisher.com/kingfisher



Reverse transcription

Reverse transcription is the synthesis, by reverse transcriptase, of complementary DNA (cDNA) using single-stranded RNA as a template. The cDNA can be used as a template for PCR amplification, cDNA library construction, RNA sequencing, and more. Selecting the right reverse transcriptase is critical to detecting low-abundance RNAs in a sample and obtaining high yields of full-length cDNA.

Find technical resources on reverse transcription at
[thermofisher.com/rteeducation](https://www.thermofisher.com/rteeducation)

Considerations for selecting the right reverse transcriptase

Sensitivity, thermostability, processivity, and inhibitor tolerance of reverse transcriptases all affect the quantity and length of cDNA synthesized.

Sensitivity

The ability of a reverse transcriptase to generate cDNA from the lowest amount of input RNA is an important attribute when working with low-copy genes or difficult sample sources where RNA may have already degraded.

Thermostability


Thermostable reverse transcriptases allow reactions to occur at higher temperatures, which help denature RNA with strong secondary structure or high GC content, for generation of longer cDNA, higher cDNA yields, and better coverage of RNA populations in the cDNA.

Processivity

Processivity is the enzyme's ability to add consecutive nucleotides without releasing the template. Highly processive reverse transcriptases allow synthesis of longer cDNA strands in a shorter reaction time, and overall better efficiency in making full-length cDNA.

Inhibitor tolerance


Compounds that have inhibitory effects on reverse transcriptases are common in RNA samples even after purification. Their sources include reagents used for RNA isolation and contaminants carried over from biological samples. Reverse transcriptases resistant to common inhibitors help minimize inconsistent or suboptimal results in cDNA-based assays.



(A) Sensitivity, (B) thermostability, (C) processivity, and (D) inhibitor tolerance of reverse transcriptases can affect the quantity and length of cDNA.

Find out more at [thermofisher.com/reverse-transcription](https://www.thermofisher.com/reverse-transcription)

Reverse transcription reagent selection guide



| Product format | Stand-alone enzyme | First-strand cDNA synthesis kit | First-strand cDNA synthesis master mix for RT-qPCR | One-step RT-PCR kit | Direct RT kit |
|---|--|--|--|--|---|
| Recommended product | Invitrogen™ SuperScript™ IV Reverse Transcriptase | Invitrogen™ SuperScript™ IV First-Strand Synthesis System | Invitrogen™ SuperScript™ IV VILO™ Master Mix | Invitrogen™ SuperScript™ IV One-Step RT-PCR System | Invitrogen™ SuperScript™ IV CellsDirect™ cDNA Synthesis Kit |
| Applications | RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq | RT-PCR, RT-qPCR, cloning, cDNA library construction, RACE, RNA-Seq | RT-qPCR | RT-PCR | RT-PCR, RT-qPCR |
| Input total | 1 pg–5 µg | 1 pg–5 µg | 0.01 pg–2.5 µg | 0.01 pg–1 µg | 1–10,000 cells |
| Optimal reaction temperature | 50–55°C | 50–55°C | 50–55°C | 50–55°C | 50–55°C |
| Reverse transcription time | 10 min | 10 min | 10 min | 10 min | 10 min |
| High cDNA yield with challenging or degraded RNA | • | • | • | • | • |



Did you know?

The standard enzyme format is incompatible for lyophilization because of the glycerol in the storage buffer. The lyo-ready (lyophilization-ready) format of SuperScript reverse transcriptases has a glycerol content below 0.1% and offers greater stability for lyophilized molecular assay kits. Learn more at thermofisher.com/lyoreadyenzymes

Find out more at thermofisher.com/superscript

Genomic DNA removal

RNA purification methods, including protocols with DNase digestion on-column, often fail to remove gDNA completely. Amplification of contaminating gDNA can cause nonspecific results. Traditional gDNA decontamination protocols with DNase I include time-consuming DNase inactivation or removal steps under conditions that can damage RNA and affect results.

Both the SuperScript IV One-Step RT-PCR System and SuperScript IV VILO Master Mix are available in a format with the novel dsDNA-specific Invitrogen™ ezDNase™ enzyme, which enables efficient, fast, and gentle (<5 min at 37°C) gDNA removal from RNA samples to help ensure high confidence in RT-PCR and RT-qPCR results.


SuperScript IV and SuperScript IV VILO Master Mix cDNA synthesis workflow with ezDNase enzyme



~27
minutes



Traditional cDNA synthesis workflow with DNase I





~105
minutes

Find out more at thermofisher.com/ssiv-onestep and thermofisher.com/4vilo

For full-length first-strand cDNA synthesis, oligo(dT) primers are recommended because of their specificity for eukaryotic mRNA, and they allow many different targets to be studied from

For target mRNA containing strong transcriptional pauses, random primers are better suited because they anneal throughout the target molecules. They are also ideal for nonpolyadenylated RNA, such as bacterial RNA.




To avoid poly(A) slippage during priming, anchored oligo(dT) primers can be used to anneal to the 5' end of the poly(A) tail of mRNA and prevent priming within the poly(A) tail. Learn more about selection of primers for reverse transcription at thermofisher.com/rteducation

Find out more at thermofisher.com/rtpprimers



PCR



The polymerase chain reaction (PCR) is a sensitive and efficient method for amplifying a single copy of a target DNA sequence to millions of copies. DNA amplification by PCR is an important step in cloning, gene expression analysis, genotyping, sequencing, and mutagenesis. PCR has a broad range of applications, including in research for infectious diseases, cancer, forensic analysis, and agricultural biotechnology.

Find technical resources on PCR at
[thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation)

Thermal cyclers

Thermal cyclers, which automate the heating and cooling cycles required to amplify DNA, play a critical role in the success of PCR. The following are things to consider when selecting a thermal cycler.

Precise temperature control

Thermal cyclers with precise temperature control enable you to quickly and accurately determine optimal annealing temperatures. Several block technologies, including gradient and Applied Biosystems™ VeriFlex™ temperature control, are available. A VeriFlex Block employs a separate heating and cooling element in each temperature zone, allowing better control and precision of temperatures. Learn more about the technology at thermofisher.com/veriflextechnology

Reliability

Thermal cyclers should be able to withstand repeated use, environmental stress, and shipping conditions. Component reliability can be tested using robotic assemblies in repeated testing of frequently used instrument components such as the heated lid, touchscreens, and temperature cycling modules. Applied Biosystems™ thermal cyclers adhere to stringent reliability criteria, which are reported at thermofisher.com/thermalcyclerreliability

Temperature accuracy

Thermal cycler temperature accuracy is a key factor in the success or failure of a PCR reaction. It is particularly important during annealing

temperature optimization, which requires both accuracy and consistency in the thermal cycler block. If the temperature set point of the instrument does not correspond to the actual temperature of the block, further temperature optimization could be required. Review a study of temperature accuracy in a number of models, available at thermofisher.com/thermalcycleraccuracy

Features

A variety of Applied Biosystems thermal cyclers are available to fit your applications and budget. Certain features may be important to you, depending on your needs. If you perform PCR optimization frequently, you will likely benefit from an instrument with a VeriFlex Block. If you would like to run optimized assays on a new or different thermal cycler, you can save re-optimization time by using a simulation mode.

If you want remote access to your instrument, you will appreciate the convenience of cloud-enabled thermal cyclers. They allow you to design and share protocols, schedule an instrument, start or stop a run, and check run status from anywhere, on any mobile device or desktop computer.


If you manage multiple thermal cyclers and users, you may benefit from a single interface for viewing all instruments at a glance and setting custom permissions by instrument, user, and method. Learn more at thermofisher.com/fleetcontrol



Helpful tip

Using the right PCR plastics for your application and instrument can improve the reliability of your PCR results. Go to thermofisher.com/findplastics to determine the right PCR plastics for you.

VeriFlex Block
temperature zones



Find out more at thermofisher.com/thermalcyclers

Select the Applied Biosystems™ thermal cycler that's right for you



Cloud-enabled



Cloud-enabled

ProFlex™ PCR System

- › Do you share the device with colleagues?
- › Do you expect your throughput needs to change?
- › Do you want to access your instrument remotely?

VeritiPro™ Thermal Cycler

- › Do you perform a lot of optimizations?
- › Do you want to access your instrument remotely?

Key benefits

Ultimate flexibility and throughput

Ultimate performance

Max sample throughput

480,000 reactions

96 reactions

Max block ramp rate

6.0°C/sec

6.0°C/sec

Temperature optimization

6-zone VeriFlex Block on 96-well system
2-zone VeriFlex Block on 3 x 32-well system

6-zone VeriFlex Block on 96-well system

Compatible with Fleet Control Software

Yes

Yes

Do you require an FDA Class 1/CE-IVD labeled thermal cycler?

Visit thermofisher.com/veritidx



Cloud-enabled



Cloud-enabled



Automation-ready

SimpliAmp™ Thermal Cycler

- › Do you need an intuitive interface?
- › Do you train new technicians often?
- › Do you want to access your instrument remotely?

Elegantly simple and precise

96 reactions

4.0°C/sec

3-zone VeriFlex Block on 96-well system

Yes

MiniAmp™ Thermal Cycler

- › Do you want an instrument with just the features needed for routine PCR?
- › Do you want to access your instrument remotely?

Routine PCR, elevated

96 reactions

3.5°C/sec

3-zone VeriFlex Block
on MiniAmp™ Plus model

Yes

Automated Thermal Cycler

- › Do you want to place your instrument on a robotic platform now or in the future?

Designed for easy robotic integration

384 reactions

3.5°C/sec

None

Yes

Find out more at thermofisher.com/thermalcyclers

PCR and qPCR plastics, seals, and accessories

Since PCR is a sensitive detection method, PCR plastics must be of high quality and free of contaminants and inhibitors, to help enable optimal performance. Regardless of the plastics format you select, proper fit and uniform heat transfer during thermal cycling are essential.

Manufacturing quality control

Applied Biosystems™ PCR and qPCR plastic consumables are manufactured in world-class facilities dedicated to the production of high-quality molecular biology-grade plastics. After manufacturing, all plastics undergo stringent quality control.

Integrity testing: Every well of every plate is visually inspected and leak tested. This thorough screening verifies every well is intact to protect all reactions.

Evaporation testing: Samples are run through PCR to test sealing performance. Well liquid volumes are analyzed post-PCR to verify seal integrity. This helps ensure that every production lot conforms to strict tolerances.

Biological testing: Our plastics are biologically tested to certify them as free of DNA, RNase, and PCR inhibitors. We offer General Purpose Laboratory Equipment (GPLE) plastics that are provided with a PCR certificate for your convenience and documentation.


Construction materials

Applied Biosystems™ MicroAmp™ optical microplates are made of polypropylene for optimal transfer of thermal energy for efficient PCR. A select medical-grade polypropylene is chosen for its exceptional biocompatibility and inert properties.

Applied Biosystems™ EnduraPlate™ microplates are constructed with a stronger polycarbonate frame to resist distortion caused by robotic grippers and to better tolerate rapid heating and cooling, while retaining thin-walled polypropylene wells for efficient heat transfer to the reaction mixture. The polycarbonate frames of the plates are available in multiple colors to help with organization and visual monitoring of assays in a high-throughput setting.



MicroAmp optical microplates.



EnduraPlate microplate for high-throughput robotic applications.

Find out more at thermofisher.com/pcrplastics

Applied Biosystems PCR and qPCR plastics are validated and tested for reliability and optimal performance. They are “Engineer Approved” for use with all Applied Biosystems thermal cyclers and real-time PCR instruments, and are available in a variety of 32-, 48-, 96-, and 384-well plates; tube strips; single tubes; caps; and seals. The table below provides a detailed comparison of each product. Easily find the PCR and qPCR plastics compatible with your instrument using the online selection tool at thermofisher.com/findplastics.

| | Small-scale experiments with a few samples | Daily experiments | Complete-workflow experiments—ideal for automation | Diagnostic procedures and automation compatible |
|----------------------------------|--|--|---|---|
| | Single tubes, strips, caps, adhesive film, and accessories | MicroAmp™ optical microplates | MicroAmp™ EnduraPlate™ optical microplates | EnduraPlate™ optical GPLE* reaction plates |
| Formats | <ul style="list-style-type: none"> Single tubes Single tubes with caps 8-strip tubes with caps 12-strip caps | <ul style="list-style-type: none"> 32-well 48-well Fast 96-well 96-well Fast 384-well | <ul style="list-style-type: none"> 96-well 96-well Fast 384-well 96-well full skirted | <ul style="list-style-type: none"> 96-well 96-well Fast 384-well |
| DNA-, RNase-, PCR inhibitor-free | Yes | Yes | Yes | Yes |
| Colors available | Clear, or mixed packs containing red, orange, blue, and green | Clear | Single-color packs (red, blue, green, yellow, or clear) and 5-plate sampler (one of each color) | Clear |
| Barcode available | No | Yes (1 or 2 sides) | Yes (3 sides) | Yes (3 sides) |
| Automation compatible | | | Yes | Yes |

* Each lot of reaction plates is certified in an ISO 13485–registered facility to be free of DNA, RNase, and PCR inhibitors. Ideal for use in diagnostic procedures.



Did you know?

Low-profile plastics, also referred to as “Fast” tubes or plates, are generally required for fast (0.1 mL) thermal blocks. Fast plastics utilize lower volumes (0.1 mL) than the standard (0.2 mL) tubes or plates. The low profile minimizes the air space above the reaction, helping reduce the effects of evaporation and enhancing thermal conductivity. Learn more about PCR and qPCR plastics at thermofisher.com/pcrplastics-education

PCR reagents

DNA polymerase is an essential component for PCR because of its key role in synthesizing new DNA strands. Because of the sensitive and specific nature of PCR, it is important to choose high-quality enzymes and reagents to produce optimal results. The following are things to consider when choosing PCR enzymes.

Specificity

Nonspecific amplification is one of the major hurdles in PCR, since it can drastically impact yield and sensitivity of target amplification. One way to help reduce nonspecific amplification is through the use of a hot-start DNA polymerase, which utilizes an antibody or chemical modification so that the polymerase becomes active only at the high temperature of the denaturation step. In addition to improving specificity, a hot-start DNA polymerase increases yield and allows convenient room temperature setup for high-throughput applications.

Thermostability

Since thermal cycling is a key feature of the conditions that enable the repetitive chain reaction of amplifying DNA, thermostability of the DNA polymerase to be used is also an important feature. Highly thermostable DNA polymerases are recommended for amplifying GC-rich or long templates that often require prolonged high-temperature reactions.

Fidelity

The fidelity, or proofreading capability, of a DNA polymerase is based on its 3' to 5' exonuclease activity, which corrects misincorporated nucleotides. This function is critical in applications such as cloning, sequencing, and site-directed mutagenesis, for accurate replication of DNA sequences.

Processivity

A DNA polymerase's processivity is defined as the number of nucleotides being incorporated in a single binding event. This property often reflects synthesis rate and speed, as well as affinity for its substrates. Therefore, highly processive DNA polymerases are beneficial to amplify challenging templates such as long, GC-rich, or inhibitor-containing DNA.

Primer annealing temperature

The primer annealing temperature of each DNA fragment to be amplified often needs optimization when designing a PCR protocol. To help simplify annealing and enable co-cycling of PCR assays, consider a DNA polymerase with a reaction buffer that allows a universal annealing temperature of 60°C for primers.



Did you know?

The residual bacterial DNA in recombinant PCR enzymes poses challenges in microbial genome analysis, such as accurately detecting bacterial strains by 16S rRNA gene sequences. To enable confidence and success in microbial PCR assays, choose PCR enzymes with controlled low levels of residual bacterial and human genomic DNA.

Find out more at

thermofisher.com/broad-range-pcr



Helpful tip

Direct PCR is a way to help simplify PCR experiments, save time, and prevent sample loss in the workflow. Direct PCR allows you to amplify target sequences directly from the samples without the need to first isolate and purify the DNA.




Find out more at

thermofisher.com/direct-pcr

Find out more at thermofisher.com/pcrenzymes

Antibody-based hot-start DNA polymerase and its activation in PCR for enhanced specificity



Helpful tip

One of the most common PCR troubleshooting issues is the presence of unwanted bands, or nonspecific amplification. To reduce nonspecific amplification:

- Optimize annealing temperature
- Check primer design
- Use hot-start PCR
- Prevent DNA cross-contamination
- Decrease template and/or primer concentration
- Optimize Mg^{2+} concentration

Choose the right PCR reagent for your research needs

A comprehensive portfolio of PCR enzymes and master mixes is available with the high performance and consistency you need. Start with the selection guide below to find the best enzyme for common PCR applications.

| DNA polymerase | Invitrogen™ Platinum™ SuperFi™ II DNA Polymerase | Invitrogen™ Platinum™ II Taq Hot-Start DNA Polymerase | Applied Biosystems™ AmpliTaq Gold™ 360 DNA Polymerase | Invitrogen™ Platinum™ Direct PCR Universal Master Mix |
|--|---|--|---|--|
| PCR type | High-fidelity PCR | Hot-start PCR | Hot-start PCR | Direct PCR |
| Capabilities | Highly accurate amplicon sequences, universal primer annealing, robust amplification of difficult targets | Universal primer annealing, fast DNA synthesis, detection of low-abundance targets | Chemical hot start | Detection of target DNA without genomic DNA purification |
| Technical specifications | | | | |
| Fidelity compared to Taq polymerase | >300x | 1x | 1x | 1x |
| Target length | Up to 20 kb* | Up to 5 kb | Up to 5 kb | Up to 8 kb |
| Hot-start modification | Antibody-mediated | Antibody-mediated | Chemical modification | Antibody-mediated |
| Speed | 15–30 sec/kb | 15 sec/kb | 60 sec/kb | 20 sec/kb |
| Universal primer annealing | Yes | Yes | | Yes |
| Inhibitor tolerance | Yes | Yes | | Yes |
| Blunt or 3'-A end | Blunt | 3'-A | 3'-A | 3'-A |
| Compatible with Applied Biosystems™ TaqMan® probes | | Yes | Yes | |
| Certified low level of bacterial gDNA | Yes | Yes | Yes | |
| Applications | | | | |
| Cloning and subcloning | • | | | |
| Site-directed mutagenesis | • | | | |
| GC-rich amplification | • | • | • | • |
| Template generation for sequencing | • | • | • | • |
| High-throughput PCR | • | • | | • |
| Long PCR (up to 20 kb) | • | | | |
| Genotyping | • | • | • | • |
| Amplification of samples with suboptimal purity | • | • | | • |
| Colony PCR | • | • | • | • |
| Multiplex PCR | • | • | • | • |
| Fast PCR | • | • | | • |

* Amplification of up to 40 kb fragment sizes is possible, but may require additional optimization of reaction conditions and primer design.

Innovations for superior PCR

PCR enzymes and reagents are continually being improved to help you get to your research destination faster. For example, the latest Platinum DNA polymerases are designed with the following key innovative features.

More robust and versatile

Advanced enzymatic engineering and methodology provide DNA polymerase with fast cycling, high tolerance of PCR inhibitors, and efficient amplification of challenging DNA like GC-rich sequences. These features help you amplify DNA targets confidently with speed and simplicity.

Find out more at thermofisher.com/platinumenzymes


Universal primer annealing

The innovative Platinum PCR buffers enable universal primer annealing at 60°C. This design allows you to co-cycle different PCR assays (instead of running them sequentially), drastically reducing tedious optimization steps and saving time.


Find out more at thermofisher.com/universalannealing

Direct gel loading

The latest Platinum DNA polymerases are available in a green buffer format that allows direct gel loading and eliminates tedious steps of dye addition, helping reduce pipetting errors. DNA migration is easily tracked with two dyes (blue and yellow) that are readily visible during electrophoresis (the lanes for 5 and 15 min in the figure to the right).



PCR



PCR primers and DNA oligos

Good design (i.e., good sequence selection) and high-quality primers are critical to your PCR reactions. In general, a length of 18–30 nucleotides for primers is optimal. The melting temperatures (T_m) of the primers should be between 65°C and 75°C, and within 5°C of each other.

OligoPerfect primer designer

Whether you are performing PCR, cloning, or capillary electrophoresis (CE) sequencing, take advantage of the benefits offered by our robust and easy-to-use Primer3-based Invitrogen™ OligoPerfect™ Designer.



Speed up—design primers for up to 50 genes at the same time



Store your data—ability to save your projects



Work smarter—recognizes .txt and .fasta file types



Order with ease—seamlessly integrates with the Invitrogen™ ordering portal

Try the OligoPerfect Designer at thermofisher.com/oligoperfect-designer

Custom DNA oligos

Invitrogen™ custom DNA oligos are synthesized on a highly automated, computer-controlled system, followed by rigorous quality control. Mass spectrometry and capillary electrophoresis are performed for short and long oligos, respectively, to help ensure the quality of the process and end products.

The appropriate synthesis scale and purification for your application depend on the nature of your downstream applications. Choose the right oligos and purification methods for your applications.



Helpful tip

If the T_m of your primer is very low, try to find a sequence with higher GC content; alternatively, the length of the primer can be extended. For more tips on primer design, go to thermofisher.com/primerdesign

| Processing option | Desalted | HPLC | PAGE |
|-------------------|---|---|--|
| Oligos | <ul style="list-style-type: none">• 25 nmol–10 μmol• 5–100 bp | <ul style="list-style-type: none">• 50 nmol–10 μmol• 7–55 bp• >85% full-length sequence | <ul style="list-style-type: none">• 50 nmol–10 μmol• 7–100 bp• >90% full-length sequence |
| Standard PCR | • | | |
| Specialty PCR | | • | • |
| Cloning | | • | • |

In addition to standard delivery, next-day delivery is also available in most regions where local oligo manufacturing exists. For more on ordering information, yield guarantees, designing tools, technical resources, protocols, and FAQs, go to thermofisher.com/oligos

Find out more at thermofisher.com/oligos



Electrophoresis

Nucleic acid electrophoresis is a common technique in molecular biology to separate, identify, quantify, and/or purify nucleic acids. Setting up electrophoresis involves a number of steps to achieve optimal separation and analysis of nucleic acid samples, such as gel preparation, ladder selection, sample visualization, and gel documentation.

Find technical resources on nucleic acid electrophoresis at
thermofisher.com/na-electrophoresis-education

Nucleic acid electrophoresis

Choosing the right tools for nucleic acid electrophoresis can significantly improve and accelerate results, enabling you to address downstream applications sooner.


Determining the appropriate gel type and gel concentration is an essential step that will help streamline the separation of nucleic acids. Learn more about convenient reagents for agarose gel electrophoresis, including hassle-free precast Invitrogen™ E-Gel™ agarose gels and pour-your-own Invitrogen™ UltraPure™ agarose reagents, in this section.

| If you need ... | Rapid results, quality control, and a safer workflow | High-quality reagents, a versatile workflow, and cost savings |
|-------------------------|---|---|
| Product | E-Gel precast agarose gels | UltraPure Agarose |
| Product format | Precast agarose cassettes | Powder |
| Protocol time (approx.) | 18 min | 120 min |
| Ready to use | Yes | No |
| Get more information at | thermofisher.com/egel | thermofisher.com/ultrapure |

E-Gel electrophoresis system workflow



Traditional DNA electrophoresis workflow



Find out more at thermofisher.com/electrophoresis

Simplify electrophoresis with E-Gel precast agarose cassettes

E-Gel precast gels

Using precast agarose gels can simplify the nucleic acid electrophoresis workflow. E-Gel precast gels are self-contained and ready for use with the agarose, electrodes, and the DNA stain packaged inside a disposable cassette. There are no gels to pour, buffers to make, staining or destaining steps to perform, or gel boxes to assemble. Just load your samples and run.

E-Gel precast gels offer excellent resolution and clarity in ≤ 18 minutes and are ideal for analyzing PCR products, restriction digests, plasmid preparations, and genotyping products. To help simplify cloning workflows, Invitrogen™ E-Gel™ CloneWell™ II gels use a double-comb design to enable recovery of purified DNA for downstream applications, without the need for additional purification kits or steps.

Find out more at thermofisher.com/egel



E-Gel DNA ladders

Accurate analysis of electrophoresis bands often depends on the DNA ladder you choose for your gel run. For optimal performance on E-Gel precast agarose gels, Invitrogen™ E-Gel™ DNA ladders are formulated with chromatography-purified DNA fragments in ready-to-use loading buffer. The chromatography-based purification method results in exceptional purity and quality of the DNA fragments, while an optimal buffer formulation reduces dye masking and helps improve ladder migration for more accurate analysis.

Find out more at thermofisher.com/egel-ladders



Helpful tip

E-Gel precast gels are available in a variety of formats for routine and high-throughput applications, with different stains (see page 36) and agarose percentages (0.8%, 1.2%, 2%, and 4%). To find the right gel for your needs, see the selection guide at thermofisher.com/egelselection

Find out more at thermofisher.com/egels

E-Gel Power Snap Electrophoresis System

To help reduce user errors and workflow time, the Invitrogen™ E-Gel™ Power Snap Electrophoresis System integrates rapid, real-time nucleic acid analysis with high-resolution image capture. The system offers:

Faster analysis—

go from sample loading to image capture in as little as 18 minutes

Simple operation—

intuitive user interface with large touchscreen and integrated operating system

Safer workflow—

minimize handling of hazardous chemicals

Find out more at thermofisher.com/powersnap



E-Gel system for high-throughput electrophoresis

High-throughput electrophoresis is often performed for high-volume analysis of PCR products, plasmid preparations, and restriction digests. Gel runs and analysis can be accelerated using the E-Gel 48 and E-Gel 96 precast gels and expandable Invitrogen™ E-Base™ Electrophoresis System. The integrated design of the E-Gel™ Mother E-Base™ and Daughter E-Base™ devices saves space and allows up to 384 samples to run at one time.

Find out more at thermofisher.com/egel-highthroughput



Electrophoresis reagents

For pouring your own agarose gels, choosing high-quality agarose, optimized DNA ladders, and improved DNA stains can help you achieve optimal electrophoresis results.

DNA stains

Detection of nucleic acid samples in gels can be improved using fluorescent dyes that are safer and/or more sensitive than ethidium bromide. The Invitrogen™ SYBR™ Safe, SYBR™ Green I, and SYBR™ Gold stains provide greater safety and/or sensitivity with lower background fluorescence than the conventional ethidium bromide stain.

Find out more at thermofisher.com/stains

UltraPure reagents for electrophoresis

Invitrogen™ UltraPure™ reagents are specifically formulated to meet your nucleic acid analysis and purification needs. UltraPure agarose and reagents are made from highly pure biochemicals for maximum reliability and superior performance.

Find out more at thermofisher.com/ultrapure

DNA ladders

Invitrogen™ DNA ladders are available in a wide variety of size ranges (10 bp to 15 kb) and formats for different applications. To create DNA ladders of superior quality, each fragment is purified individually using proprietary chromatography-based technology. Our DNA ladders are stable during prolonged storage at room temperature and after multiple freeze-thaw cycles.

Find out more at thermofisher.com/ladders



Did you know?

Chromatographically purified nucleic acid fragments are considered the gold standard for ladders, since the technology provides higher control over quality, banding pattern, intensity, and quantity for ladder composition.

Learn more at thermofisher.com/na-electrophoresis-education

Fluorescent nucleic acid gel stains

| | Standard detection | Safer detection | Enhanced detection | Ultimate detection |
|--|----------------------------|------------------|---------------------------|-------------------------|
| | UltraPure ethidium bromide | SYBR Safe stain | SYBR Green I stain | SYBR Gold stain |
| Sensitivity (dsDNA) | Sensitive (1 ng) | Sensitive (1 ng) | Highly sensitive (>60 pg) | Ultrasensitive (>25 pg) |
| Less hazardous and more environmentally friendly | | • | | |
| Improved cloning efficiency | | • | • | • |



Cloning

Molecular cloning involves recombinant DNA technologies that insert a DNA sequence of interest into a vector to generate a large number of copies. Traditionally, cloning has been carried out with restriction enzymes and a DNA ligase to form a new vector capable of expressing the gene of interest. In the case of gene synthesis, researchers can obtain their desired DNA directly in a specified vector with just sequence information. Other cloning methods, such as PCR cloning, Invitrogen™ TOPO™ cloning, ligation-independent cloning, and gene assembly are commonplace, exploiting unique characteristics of other DNA-modifying enzymes.

Find technical resources on molecular cloning at
thermofisher.com/cloningeducation

Cloning and gene synthesis

From restriction enzymes to gene synthesis, a large portfolio of tools and resources is available to help you obtain high-quality cloned DNA for your next discovery.

| Method | Thermo Scientific™ FastDigest™ restriction enzymes | Invitrogen™ TOPO™ cloning | Invitrogen™ Gateway™ cloning | Invitrogen™ GeneArt™ seamless cloning and GeneArt™ Gibson Assembly® cloning kits | Invitrogen™ GeneArt™ Type IIIs assembly | Invitrogen™ GeneArt™ Strings™ DNA Fragments | Invitrogen™ GeneArt™ Gene Synthesis |
|---|---|---|--|--|--|---|--|
| Key benefits/description | <ul style="list-style-type: none"> • Familiarity, flexibility, convenience, time savings • Universal protocol with 15 min digestion in one buffer | <ul style="list-style-type: none"> • >95% efficiency, 5 min PCR cloning • Compatible with many other cloning systems | <ul style="list-style-type: none"> • High-throughput and high-efficiency shuttling among multiple expression vectors | <ul style="list-style-type: none"> • Seamless multi-fragment assembly by homologous recombination • Directional cloning of up to 15 fragments • Up to 95% efficiency and 15 min cloning | <ul style="list-style-type: none"> • One-tube seamless multi-fragment assembly by simultaneous restriction digestion and ligation • Directional cloning of up to 8 fragments, for up to 20 kb total • Efficient for repetitive and very small sequences | <ul style="list-style-type: none"> • Synthesized DNA fragments ready to clone via the method of your choice • No starting DNA required • Pool sequence-verified | <ul style="list-style-type: none"> • Custom-cloned genes in your choice of vector • Sequence-verified • Can be optimized for a specific host for maximal protein expression |
| Technology basics | <ul style="list-style-type: none"> • Restriction digestion and ligation | <ul style="list-style-type: none"> • Topoisomerase-based, ligase-free cloning | <ul style="list-style-type: none"> • Single-step, directional, and site-specific DNA recombination • Restriction enzyme- and ligase-free | <ul style="list-style-type: none"> • End-terminal homology recombination using overlapping sequences • Transformation-associated recombination (TAR) in <i>Saccharomyces cerevisiae</i> | <ul style="list-style-type: none"> • Type IIIs restriction and ligation in a single reaction | <ul style="list-style-type: none"> • Linear dsDNA assembled from pooled synthetic oligonucleotides • 150–3,000 bp, also available in library format with randomized bases | <ul style="list-style-type: none"> • DNA of interest cloned in vector • 100% sequence-verified with quality assurance documentation |
| Needs DNA source material (gene in plasmid, library, etc.) | • | • | • | • | • | | |
| Use your own vector | • | | * | • | • | • | • |

* Vector needs to be converted with Invitrogen™ Gateway™ Vector Conversion System with One Shot™ *ccdB* Survival™ 2 T1^R Competent Cells.

Find out more at thermofisher.com/cloning

Restriction enzyme cloning


Found naturally in bacteria, restriction enzymes recognize and cleave specific DNA sequences, resulting in sticky ends (5' or 3' protruding ends) or blunt ends, enabling DNA inserts to be cloned into vectors with compatible ends. Star activity, buffer compatibility, and varying protocols for complete digestion are some common hurdles in restriction digestion.

FastDigest restriction enzymes

To simplify cloning, we offer FastDigest enzymes—an advanced line of restriction enzymes that share buffer compatibility with downstream modifying enzymes. Its benefits include:

- Complete digestion in 5–15 min
- Double and multiple digestions in a universal buffer for any combination of enzymes
- No sequential digestions and buffer changes
- 176 unique specificities
- Direct loading of reaction mixture on gels

Find out more at thermofisher.com/fastdigest




Type IIs restriction enzymes

A specific group of restriction enzymes called Type IIs endonucleases cleave DNA outside of their recognition sequences. In combination with DNA ligase, Type IIs restriction enzymes are utilized to drive the insertion of one or several DNA fragments into a recipient vector without the inclusion of residual restriction enzyme sites and other unwanted DNA sequences at fragment junctions (scarless cloning).

Find FastDigest Type IIs enzymes at
thermofisher.com/fastdigesttypeiis

For GeneArt Type IIs Assembly Kits, go to
thermofisher.com/typeiis




PCR cloning

PCR cloning is a method in which double-stranded DNA fragments amplified by PCR are ligated into a vector. With PCR amplification, this cloning technique requires much less starting material for the insert sequence and allows introduction of new restriction and/or recombination sites to the 5' end of the inserts.

TOPO cloning

TOPO PCR cloning technology was developed to help improve cloning efficiency, simplify protocol setup, and accommodate a wide range of PCR insert sizes. TOPO cloning vectors are linearized by the activity of topoisomerase I (which also has a ligase function) that is covalently bound to the 3' phosphate on each end (see figure below). This system enables the vectors to readily be joined to PCR inserts with compatible ends (with up to 95% efficiency), without the need for additional ligation steps, in 5 minutes.

Find out more at thermofisher.com/topo




Did you know?

The Invitrogen™ TOPO™ XL-2 Complete PCR Cloning Kit provides all the necessary elements for highly efficient cloning of extra-long PCR products from 1–13 kb. thermofisher.com/topoxl2

Gateway cloning

To shuttle a PCR insert among vectors, the Gateway cloning system offers site-specific, recombinase-based cloning. It maintains the insert's proper orientation and reading frame during shuttling using the Gateway vectors. Once a gene is cloned into an entry clone, you can then move the DNA fragment into one or more destination vectors simultaneously.

Find out more at thermofisher.com/gateway




Seamless cloning and GeneArt Gibson Assembly cloning kits

To assemble multiple PCR fragments by end-terminal homologous recombination, several seamless cloning technologies are available for scarless and directional cloning into any vector. GeneArt seamless cloning kits offer the option of building constructs using *E. coli* and *Saccharomyces cerevisiae*.

Invitrogen™ GeneArt™ Gibson Assembly® kits allow for the simultaneous assembly of up to 15 very large DNA fragments to create precise constructs with no additional sequences, in highly efficient reactions. This cloning method circumvents the need for multiple rounds of restriction enzyme analysis and digestion, DNA end repair, dephosphorylation, ligation, enzyme inactivation, and cleanup, and is a powerful tool in synthetic biology.

GeneArt Gibson Assembly kits offer these benefits:

- Assembly of up to 15 fragments to build seamless clones
- Cloning efficiencies up to >95%
- Choice of complete kits with competent cells or master mixes



Did you know?

The Gibson Assembly method has been referenced in thousands of peer-reviewed publications and is a powerful method that can be used to seamlessly construct synthetic and natural genes, genetic pathways, and entire genomes.¹

Find out more at thermofisher.com/seamless

1. Enzymatic assembly of DNA molecules up to several hundred kilobases. Gibson DG et al. (2009) *Nat Methods* 6(5):343-5.

Cloning with synthetic DNA

If you lack the time to generate and clone insert DNA, including optimization and troubleshooting, our synthetic DNA fragments and cloning service might be right for you. GeneArt Strings DNA Fragments and GeneArt Gene Synthesis offer genes analogous to optimized, error-free PCR products.

GeneArt Strings DNA Fragments

A time-saving alternative to PCR, GeneArt Strings DNA Fragments are available in lengths up to 3 kb and are compatible with any downstream cloning method of choice, providing:

- Synthetic, ready-to-use DNA fragments
- DNA with your specified ends to facilitate the cloning method of choice
- No starting DNA required
- Free optimization of gene with Invitrogen™ GeneArt™ GeneOptimizer™ software for maximum protein expression
- Option of Strings DNA Libraries with mixed, randomized nucleotides using full IUPAC code


Find out more at thermofisher.com/strings

GeneArt Gene Synthesis

A reliable and cost-effective method for obtaining customized DNA constructs with 100% sequence accuracy, GeneArt Gene Synthesis offers:

- Synthetic, ready-to-transfect genes
- Cloning into several available vectors (custom options available)
- 100% sequence-verified and ready for downstream applications
- No starting DNA required
- Free optimization of gene with GeneOptimizer software for maximum protein expression


Find out more at thermofisher.com/genesynthesis



Transformation

Once the DNA fragment is cloned into a vector, transformation into bacteria is performed to enable propagation of sufficient quantities of the cloned DNA for downstream experiments. Selection of competent cells for transformation depends upon the transformation methods, strain genotypes, plasmid characteristics, and desired applications. Visit thermofisher.com/compcells-education for technical resources on competent cells.

Choosing Invitrogen™ competent cells based on the application



Find out more at thermofisher.com/compcells

Transformation (cont.)

Medium- and high-throughput transformation

Performing bacterial transformations one by one can be very time-consuming and create a bottleneck in your experimental workflow. There are times when medium- and high-throughput transformation options are desired. Invitrogen™ MultiShot™ chemically competent cells provide three flexible product formats to meet your throughput needs.

Find out more at thermofisher.com/multishot



StripWell format

- Medium-throughput option
- Twelve 8-tube strips
- Suitable for 1–96 transformations
- Five *E. coli* strains available

FlexPlate format

- High-throughput option
- 96-well plate separates into 12 x 8-well segments
- Manual and automated platform transformations
- Six *E. coli* strains available



96-well plate

- Highest-throughput option
- Five 96-well plates
- Available with the TOP10 strain
- Stable replication of high copy number plasmids



Did you know?

Invitrogen competent cells can be provided in custom configurations per your request. Large and custom volumes as well as multiple formats are at your fingertips. Simply email us at customorders@thermofisher.com



Resources


We have compiled educational resources, mobile apps, frequently asked questions, and other information to help you reach new heights in your research.



Educational resources

Suitable for new and experienced molecular biologists alike, our free online education resources are designed to help you review the basics, build your expertise, or discover our latest innovative technologies. Explore our technical resources in the following areas of molecular biology.


Reverse transcription




PCR enzymes, plastics, and thermal cyclers



Nucleic acid electrophoresis



Restriction enzymes, molecular cloning, and competent cells



Resources

Webinars: Watch live and recorded webinars for in-depth understanding of molecular and synthetic biology techniques and tools to help elevate your research.
thermofisher.com/mbwebinars

Videos: Experience entertaining and visual learning with our educational videos on molecular biology techniques, how-tos, tips and tricks, and more.
thermofisher.com/mbvideos

Application notes: Read white papers and application notes from our R&D scientists on our product innovations.
thermofisher.com/mbliterature

Online tools: Use our interactive online tools for PCR annealing temperature, restriction enzyme information, product selection, and more.
thermofisher.com/mbtools

Find out more at thermofisher.com/mbschool

Mobile apps



DailyCalcs—science calculator

The DailyCalcs app turns your phone into a science calculator to help simplify everyday tasks in the lab. The app features eight calculators: molarity, dilution, formula weight, transfection, unit conversions, culture vessel data, media conversions, and specific productivity.



Instrument Connect—remote monitoring

Instrument Connect allows you to view instrument status, monitor or schedule a run, and more on any cloud-enabled instrument, including the ProFlex, SimpliAmp, and MiniAmp PCR instruments.



PCR Quest—match-3 lab game


Test your PCR knowledge with our lab game—PCR Quest—where you travel from lab to lab crushing the world's toughest diseases. Download at thermofisher.com/pcrquest

Custom and OEM solutions

As a leading supplier of molecular biology reagents, plastics, and instruments, we offer customizable manufacturing solutions used by companies in developing next-generation molecular assays. Regardless of where you are in your assay development, we have off-the-shelf or custom solutions to help you achieve your goals. Partner with an experienced supplier that understands both raw materials and new technologies—a market leader with a dedicated diagnostics partnering business that brings value beyond products.

What do our OEM solutions mean to you?

- Customization of products and services
- Consultation, partnership, and expertise
- Negotiated business terms
- Warranties and indemnification
- Commercial-use rights and obligations
- Risk and liability management



Find out more at thermofisher.com/oemmolecular

Find out more at thermofisher.com/mobileapps

Frequently asked questions

Below are some common questions and answers to help you start or troubleshoot molecular biology experiments.

Sample preparation

Which kit should I use to isolate nucleic acids from my sample?

Choosing the right product is fundamental to ensuring proper lysis of cells and tissue, as well as sufficient yield and quality of isolated nucleic acids. Look to our selection guides (see pages 9–11) to help you decide according to nucleic acid type, sample source, experimental throughput, and format as well as downstream applications.

What are the key steps to preventing RNA degradation?

The basic lab precautions listed below can help minimize RNA degradation and avoid experimental inconsistency and failure.

- Use nuclease-free pipette tips and tubes
- Use nuclease-free water and reagents
- Regularly decontaminate work surfaces
- Properly stabilize RNA sources before storage

For more tips and troubleshooting advice on sample prep, visit [thermofisher.com/rnabasics](https://www.thermofisher.com/rnabasics) and [thermofisher.com/napsupport](https://www.thermofisher.com/napsupport)

Reverse transcription

How do I improve the efficiency of cDNA synthesis when working with challenging samples (e.g., low-abundance, degraded, inhibitor-containing, or GC-rich RNA)?

When working with challenging RNA samples, select a reverse transcriptase that is highly sensitive, processive, thermostable, and resistant to common inhibitors, to help you obtain the highest cDNA yield (see page 15).

What are the benefits of using random primers, oligo(dT) primers, gene-specific primers, or oligo(dT)/random mixed primers in reverse transcription?

- Random primers are good to use with degraded RNA, RNA with high secondary structure, nonpolyadenylated RNA, or prokaryotic RNA.
- Oligo(dT) primers are an optimal choice for synthesis of full-length cDNA from eukaryotic mRNA. Applications include cDNA cloning, cDNA library construction, and 3' rapid amplification of cDNA ends (3' RACE).
- Gene-specific primers are designed based on known sequences of the target RNA. These primers offer the most specific priming and are commonly used in one-step RT-PCR.
- A mixture of oligo(dT) and random primers is often used in two-step RT-PCR to achieve the benefits of each primer type (see page 18).

For more tips and troubleshooting advice on reverse transcription, visit [thermofisher.com/rteducation](https://www.thermofisher.com/rteducation) and [thermofisher.com/rtsupport](https://www.thermofisher.com/rtsupport)

PCR amplification

How can I optimize primer annealing for PCR?

Traditionally, gradient thermal cyclers have been used to simultaneously assess a number of temperatures around the theoretical annealing point. Compared to gradient thermal cyclers, instruments with the VeriFlex technology allow more precise temperature control for faster optimization of primer annealing (see page 21).

Tedious optimization steps may be circumvented using the novel Platinum DNA polymerases. Their innovative buffers enable specific annealing at 60°C for most primers when they are designed following general primer design rules (see pages 28–29).

Frequently asked questions (cont.)

What do I need to run fast PCR?

PCR amplicons shorter than 1 kb can be amplified in as little as 40 minutes using “fast” enzymes (high processivity; see page 28), “fast” plastics (low profile and ultra-thin walls; see page 25), and “fast” thermal cyclers (fast ramp rate; see pages 22–23).

How can I prevent sample evaporation during PCR?

Proper sealing of your reactions will help prevent evaporation during PCR.

- When using adhesive film to seal a plate, be sure to properly align the seal to cover all wells and press firmly along all edges of the plate using an applicator tool.
- When sealing a plate using cap strips, ensure that the cap strips are compatible with the plate and thermal cycler being used. Be sure to align cap strips with each well of the plate and place firmly across the plate for a secure fit.
- Use the applicator tool (Cat. No. 4333183 or 4330015) or other comparable sealing tools as needed.

For more tips and troubleshooting advice on PCR, visit [thermofisher.com/pcreducation](https://www.thermofisher.com/pcreducation) and [thermofisher.com/pcrsupport](https://www.thermofisher.com/pcrsupport)

Nucleic acid electrophoresis

Why is it important to choose the right ladder when using E-Gel precast agarose gels?

Accurate analysis of electrophoresis bands often depends on the DNA ladder chosen for your gel run. E-Gel DNA ladders are formulated with ready-to-use buffers unique for E-Gel precast agarose gels, and DNA standards designed for optimal separation (see page 34).

Are there safer alternatives to ethidium bromide for staining nucleic acids in gel electrophoresis?

SYBR Safe DNA gel stain is a safer alternative to ethidium bromide and is commonly used in gel electrophoresis. SYBR Safe DNA gel stain is not classified as hazardous waste or as a pollutant under US federal regulations (see page 36).

For more tips and troubleshooting advice on nucleic acid electrophoresis, visit [thermofisher.com/na-electrophoresis-education](https://www.thermofisher.com/na-electrophoresis-education) and [thermofisher.com/na-electrophoresis-support](https://www.thermofisher.com/na-electrophoresis-support)

Cloning

Do you have a buffer compatibility chart for restriction enzymes?

All FastDigest restriction enzymes are 100% active in one universal FastDigest buffer (see page 40). Hence, there is no buffer compatibility chart for FastDigest restriction enzymes.

What is the main difference between GeneArt Strings DNA Fragments and GeneArt Gene Synthesis?

GeneArt Strings DNA Fragments are custom-made, uncloned, double-stranded linear DNA fragments. GeneArt Gene Synthesis is a service offered for chemical synthesis, cloning, and sequence verification of genetic sequences (see page 44).

What are some key considerations for choosing competent cells for my cloning applications?

Genotype, transformation efficiency, growth rate, and throughput format are important factors in choosing competent cells for cloning. The genotype of a cell strain may determine growth conditions and suitability for transformation with specific DNA types (see page 45).

For more tips and troubleshooting advice on cloning, visit [thermofisher.com/cloningeducation](https://www.thermofisher.com/cloningeducation) and [thermofisher.com/cloningsupport](https://www.thermofisher.com/cloningsupport)

Ordering information

| | Quantity | Cat. No. |
|--|---------------|-----------|
| Nucleic acid isolation | | |
| PureLink Quick Plasmid Miniprep Kit | 50 preps | K210010 |
| PureLink HiPure Plasmid Filter Midiprep Kit | 25 preps | K210014 |
| PureLink HiPure Plasmid Maxiprep Kit | 10 preps | K210006 |
| PureLink <i>Pro</i> Quick96 Plasmid Purification Kit | 4 x 96 preps | K211004A |
| PureLink Quick Gel Extraction Kit | 50 preps | K210012 |
| TRIzol Plus RNA Purification Kit | 50 preps | 12183555 |
| PureLink RNA Mini Kit | 10 preps | 12183020 |
| PureLink Genomic DNA Mini Kit | 10 preps | K182000 |
| PureLink <i>Pro</i> 96 Genomic DNA Mini Kit | 4 x 96 preps | K182104A |
| PureLink <i>Pro</i> 96 Viral RNA/DNA Purification Kit | 4 plates | 12280096A |
| PureLink Viral RNA/DNA Mini Kit | 50 preps | 12280050 |
| PureLink Genomic Plant DNA Purification Kit | 50 preps | K183001 |
| MagMAX DNA Multi-Sample Ultra Kit | 500 preps | A25597 |
| KingFisher Flex Purification System with 96 Deep-Well Head | 1 system | 5400630 |
| KingFisher Duo Prime Purification System | 1 system | 5400110 |
| PureLink PCR Purification Kit | 50 preps | K310001 |
| PureLink Quick Gel Extraction and PCR Purification Kit | 50 preps | K220001 |
| PureLink Quick Gel Extraction Kit | 50 preps | K210012 |
| Dynabeads M-270 Streptavidin | 2 mL | 65305 |
| Dynabeads MyOne Streptavidin C1 | 2 mL | 65001 |
| Reverse transcription | | |
| SuperScript IV Reverse Transcriptase | 2,000 units | 18090010 |
| | 10,000 units | 18090050 |
| SuperScript IV First-Strand Synthesis System | 50 reactions | 18091050 |
| | 200 reactions | 18091200 |
| SuperScript IV VILO Master Mix | 50 reactions | 11756050 |
| | 500 reactions | 11756500 |
| SuperScript IV VILO Master Mix with ezDNase Enzyme | 50 reactions | 11766050 |
| | 500 reactions | 11766500 |
| SuperScript IV One-Step RT-PCR System | 25 reactions | 12594025 |
| | 100 reactions | 12594100 |
| SuperScript IV One-Step RT-PCR System with ezDNase Enzyme | 25 reactions | 12595025 |
| | 100 reactions | 12595100 |

| | Quantity | Cat. No. |
|---|--------------------------|----------|
| Reverse transcription (cont.) | | |
| SuperScript IV CellsDirect cDNA Synthesis Kit | 50 reactions | 11750150 |
| | 500 reactions | 11750350 |
| SuperScript IV CellsDirect Lysis Reagents | 500 reactions | 11750550 |
| RNaseOUT Recombinant Ribonuclease Inhibitor | 5,000 units | 10777019 |
| Ribonuclease H | 30 units | 18021014 |
| Random Hexamers (50 µM) | 5 nmol | N8080127 |
| Random Primers | 9 A ₂₆₀ units | 48190011 |
| Oligo(dT) ₁₂₋₁₈ Primer | 25 µg | 18418012 |
| Oligo(dT) ₂₀ Primer | 15 µg | 18418020 |
| DNase I, Amplification Grade | 100 units | 18068015 |
| PCR | | |
| DNA Oligo, Desalted, Dry | 25 nmol | A15612 |
| DNA Oligo, Desalted, Dry, next-day (ordered before 1 PM Eastern Time) | 25 nmol | A15613 |
| DNA Oligo, Desalted, Liquid | 25 nmol | A15611 |
| DNA Oligo, Desalted, Dry | 50 nmol | A15610 |
| DNA Oligo, Desalted, Liquid | 50 nmol | A15609 |
| DNA Oligo, Cartridge, Dry | 50 nmol | A15614 |
| DNA Oligo, Cartridge, Liquid | 50 nmol | A15608 |
| DNA Oligo, HPLC, Dry | 50 nmol | A15607 |
| DNA Oligo, HPLC, Liquid | 50 nmol | A15606 |
| DNA Oligo, PAGE, Dry | 50 nmol | A15605 |
| DNA Oligo, PAGE, Liquid | 50 nmol | A15604 |
| Platinum II <i>Taq</i> Hot-Start DNA Polymerase | 100 reactions | 14966001 |
| | 500 reactions | 14966005 |
| Platinum II Hot-Start PCR Master Mix (2X) | 50 reactions | 14000012 |
| | 200 reactions | 14000013 |
| Platinum II Hot-Start Green PCR Master Mix (2X) | 50 reactions | 14001012 |
| | 200 reactions | 14001013 |
| AmpliTaq Gold 360 DNA Polymerase | 100 units | 4398813 |
| | 250 units | 4398823 |
| AmpliTaq Gold 360 Master Mix | 1 mL | 4398876 |
| | 5 mL | 4398881 |

Ordering information (cont.)

| | Quantity | Cat. No. |
|---|---------------|-----------|
| PCR (cont.) | | |
| Platinum SuperFi II DNA Polymerase | 100 units | 12361010 |
| | 500 units | 12361050 |
| Platinum SuperFi II PCR Master Mix | 100 reactions | 12368010 |
| | 500 reactions | 12368050 |
| Platinum SuperFi II Green PCR Master Mix | 100 reactions | 12369010 |
| | 500 reactions | 12369050 |
| Platinum Direct PCR Universal Master Mix | 100 reactions | A44647100 |
| | 500 reactions | A44647500 |
| dNTP Set (100 mM) | 4 x 250 µL | 10297018 |
| | 8 x 1.25 mL | 10297117 |
| ProFlex 3 x 32-Well PCR System | 1 instrument | 4484073 |
| ProFlex 96-Well PCR System | 1 instrument | 4484075 |
| SimpliAmp Thermal Cycler | 1 instrument | A24811 |
| VeritiPro 96-Well Thermal Cycler | 1 instrument | A48141 |
| MiniAmp Plus Thermal Cycler | 1 instrument | A37835 |
| MiniAmp Thermal Cycler | 1 instrument | A37834 |
| Automated Thermal Cycler, 96-well | 1 instrument | A31486 |
| MicroAmp EnduraPlate Optical 96-Well Fast Multicolor Reaction Plates with Barcode | 5 plates | 4483493 |
| MicroAmp Optical Adhesive Film | 100 covers | 4311971 |
| MicroAmp Optical 96-Well Reaction Plate | 10 plates | N8010560 |
| MicroAmp Optical 8-Cap Strips | 300 strips | 4323032 |
| MicroAmp Fast Optical 96-Well Reaction Plate, 0.1 mL | 10 plates | 4346907 |
| MicroAmp Fast Reaction Tube with Cap, 0.1 mL | 1,000 tubes | 4358297 |
| MicroAmp EnduraPlate Optical 384-Well Multicolor Reaction Plates with Barcode | 5 plates | 4483316 |
| MicroAmp EnduraPlate Optical 96-Well Clear Reaction Plates with Barcode | 20 plates | 4483354 |
| MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate | 20 plates | A32811 |
| MicroAmp 8-Tube Strip with Attached Domed Caps, 0.2 mL | 125 strips | A30589 |
| MicroAmp EnduraPlate Optical 96-Well Full-Skirted Plates with Barcode, clear | 50 plates | A31728 |

| | Quantity | Cat. No. |
|---|------------------|----------|
| Nucleic acid separation and analysis | | |
| UltraPure Ethidium Bromide, 10 mg/mL | 10 mL | 15585011 |
| SYBR Safe DNA Gel Stain | 400 µL | S33102 |
| SYBR Gold Nucleic Acid Gel Stain | 500 µL | S11494 |
| UltraPure DNase/RNase-Free Distilled Water | 500 mL | 10977015 |
| UltraPure Agarose | 100 g | 16500100 |
| TrackIt 100 bp Plus DNA Ladder | 100 applications | 10488058 |
| UltraPure TAE Buffer, 10X | 4 L | 15558026 |
| E-Gel Agarose Gels with SYBR Safe stain | 10 gels | A42135 |
| E-Gel Double Comb Agarose Gels with SYBR Safe stain | 10 gels | A42348 |
| E-Gel EX Double Comb Agarose | 10 gels | A42346 |
| E-Gel CloneWell II Agarose Gels with SYBR Safe DNA Gel Stain, 0.8% | 18 gels | G661818 |
| E-Gel Agarose Gels with SYBR Safe DNA Gel Stain, 2% | 18 gels | G521802 |
| E-Gel EX Agarose Gels, 2% with SYBR Gold DNA stain | 20 gels | G402002 |
| E-Gel 1 Kb Plus DNA Ladder | 100 applications | 10488090 |
| E-Gel Sample Loading Buffer, 1X | 4 x 1.25 mL | 10482055 |
| E-Gel Power Snap Electrophoresis System Starter Kit, EX 2% | 1 kit | G8342ST |
| E-Gel Power Snap Electrophoresis Device Starter Kit, CloneWell II 0.8% with SYBR Safe gel stain | 1 kit | G8168ST |
| E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 2% | 8 gels | G820802 |
| E-Gel 96 Agarose Gels with SYBR Safe DNA Gel Stain, 2% | 8 gels | G720802 |
| E-Gel 48 Agarose Gels with SYBR Safe DNA Gel Stain, 4% | 8 gels | G820804 |
| Cloning and gene synthesis | | |
| FastDigest BamHI | 800 reactions | FD0054 |
| | 2,500 reactions | FD0055 |
| FastDigest BclI | 20 reactions | FD1253 |
| | 50 reactions | FD1254 |
| FastDigest BshTI | 20 reactions | FD1464 |
| | 50 reactions | FD1703 |
| FastDigest DpnI | 100 reactions | FD1704 |
| | 800 reactions | FD0274 |
| FastDigest EcoRI | 2,500 reactions | FD0275 |
| | 300 reactions | FD0524 |

| | Quantity | Cat. No. |
|---|-----------------|----------|
| Cloning and gene synthesis (cont.) | | |
| FastDigest NotI | 20 reactions | FD0593 |
| | 50 reactions | FD0594 |
| | 150 reactions | FD0595 |
| | 250 reactions | FD0596 |
| FastDigest Sall | 200 reactions | FD0644 |
| FastDigest XbaI | 300 reactions | FD0684 |
| | 750 reactions | FD0685 |
| FastDigest XhoI | 400 reactions | FD0694 |
| | 1,200 reactions | FD0695 |
| FastDigest Esp3I (BsmBI) (IIs class) | 20 reactions | FD0454 |
| FastDigest BpiI (BbsI) (IIs class) | 20 reactions | FD1014 |
| FastDigest Eco31I (BsaI) (IIs class) | 50 reactions | FD0293 |
| | 100 reactions | FD0294 |
| TOPO TA Cloning Kit for Subcloning, without competent cells | 25 reactions | 450641 |
| Zero Blunt TOPO PCR Cloning Kit, without competent cells | 25 reactions | 450245 |
| pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i> | 20 reactions | K240020 |
| pcDNA 6.2/V5-PL-DEST Mammalian Expression Vector | 6 µg | 12537162 |
| One Shot TOP10 Chemically Competent <i>E. coli</i> | 20 reactions | C404003 |
| One Shot Stbl3 Chemically Competent <i>E. coli</i> | 20 x 50 µL | C737303 |
| MAX Efficiency DH5α Competent Cells | 200 µL | 18258012 |
| ElectroMAX DH10B Cells | 100 µL | 18290015 |
| MAX Efficiency Stbl2 Competent Cells | 5 x 200 µL | 10268019 |
| MultiShot TOP10 Chemically Competent <i>E. coli</i> | 5 plates | C40005 |
| MultiShot StripWell TOP10 Chemically Competent <i>E. coli</i> | 1 rack | C409601 |
| MultiShot StripWell BL21 Star (DE3) Chemically Competent <i>E. coli</i> | 1 rack | C609601 |
| MultiShot FlexPlate TOP10 Chemically Competent <i>E. coli</i> | 1 plate | C4081201 |
| MultiShot FlexPlate DH5α T1 ^R Chemically Competent <i>E. coli</i> | 1 plate | C4481201 |
| MultiShot FlexPlate Stbl3 Chemically Competent <i>E. coli</i> | 1 plate | C7381201 |

| | Quantity | Cat. No. |
|---|--------------|----------|
| Cloning and gene synthesis (cont.) | | |
| GeneArt Gibson Assembly HiFi Cloning Kit, chemically competent cells | 10 reactions | A46624 |
| GeneArt Gibson Assembly EX Cloning Kit, chemically competent cells | 10 reactions | A46633 |
| GeneArt Seamless PLUS Cloning and Assembly Kit | 20 reactions | A14603 |
| GeneArt Type IIs Assembly Kit, AarI | 10 reactions | A15916 |
| GeneArt Type IIs Assembly Kit, BsaI | 10 reactions | A15917 |
| GeneArt Type IIs Assembly Kit, BbsI | 10 reactions | A15918 |
| GeneArt High-Order Genetic Assembly System | 10 reactions | A13285 |
| Gateway BP Clonase II Enzyme Mix | 20 reactions | 11789020 |
| Gateway LR Clonase II Enzyme Mix | 20 reactions | 11791020 |
| MultiSite Gateway Pro Plus | 20 reactions | 12537100 |
| LR Clonase II Plus Enzyme | 20 reactions | 12538120 |
| Gateway Vector Conversion System with One Shot <i>ccdB</i> Survival Cells | 1 kit | 11828029 |
| PCR Cloning System with Gateway Technology with pDONR 221 and OmniMAX 2 Competent Cells | 20 reactions | 12535029 |
| PCR Cloning System with Gateway Technology with pDONR/Zeo and OmniMAX 2 Competent Cells | 20 reactions | 12535037 |
| Gateway pDONR 221 Vector | 6 µg | 12536017 |
| pENTR/D-TOPO Cloning Kit, with One Shot TOP10 Chemically Competent <i>E. coli</i> | 20 reactions | K240020 |
| pCR 8/GW/TOPO TA Cloning Kit with One Shot TOP10 <i>E. coli</i> | 20 reactions | K250020 |

GeneArt Gene Synthesis thermofisher.com/genesynthesis

GeneArt Strings DNA Fragments thermofisher.com/strings

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

In Canada

Order online: fishersci.ca

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437



For Research Use Only. Not for use in diagnostic procedures. © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Ensembl is a trademark of Genome Research Limited. TaqMan is a registered trademark and AmpliTaq Gold is a trademark of Roche Molecular Systems, Inc., used under permission and license. TRIzol and DNAzol are trademarks of Molecular Research Center, Inc. Gibson Assembly is a registered trademark of SGI-DNA, Inc. NEB is a trademark of New England Biolabs Inc. **BN20200633 0920**